

Focus On ALIGNER SYSTEMS

3M Oral Care
3m.com/clarityaligners
(800) 276-8789





Align Technology Inc
invisalign.com/partner
(888) 822-5446

Allesee Orthodontic Appliances
aaaaccess.com
(800) 262-5221

Dentsply Sirona
suresmile.com
(888) 672-6387



Information in this guide based on data submitted by product manufacturers. Orthodontic Products strives for accuracy in all data but cannot be held responsible for claims made by manufacturers.

				
Product name	3M Clarity Aligners	Invisalign Clear Aligners	Simpli5	SureSmile Aligner
What materials is the product made of?	Plastic	Invisalign SmartTrack material—a proprietary multi-layered polymer designed specifically for use with the Invisalign system.	Proprietary	Essix plastics
Made in lab or in-house	In-house	In-house	Lab	Both
Software interface	Access the HIPAA compliant 3M Oral Care portal at any time and from anywhere to easily refine and update treatment plans. Oral Care portal is an open platform allowing you the flexibility to use the intraoral scanner of your choice. Choose from a variety of esthetic treatment options from 3M to ensure optimal care for your patient.	ClinCheck Pro with 3D Controls	N/A	SureSmile software is cloud based for easy access and delivers a complete, clinician controlled clear aligner treatment solution for both full service and in-office fabrication.
Orthodontist's cost per case	Contact company	Varies by product	Contact company. Discounted pricing on single arch cases.	Varies with product and annual volume
Special features	The 3M Clarity Aligner system delivers custom esthetic treatment that showcases orthodontists' expertise. These removable, nearly invisible, clear aligners deliver control, flexibility and choice. In conjunction with Clarity Aligners, 3M offers the HIPAA-compliant 3M Oral Care portal, a web-based portal allowing access to treatment plans anytime and anywhere.	Weekly aligner changes may reduce treatment time up to 50% compared to 2-week changes. Precision cuts for use of elastics—designed to eliminate need to manually cut aligners. Depending on package, additional aligners may be available at no additional charge for up to 5 years maximum per case.	Simple and reliable, Simpli5 has been a dependable option for Express aligner treatment for anyone looking to perfect a beautiful smile. No approver required.	Clinically proven for over 15 years, the SureSmile software is able to link CBCT/Ceph/Facial photos with a 3D model to ensure predictable and biologically achievable results with optimal smile design. Treatment plans are digitally designed with ABO graded finish. Automated attachment placement and staging are fully adjustable by the clinician.






DynaFlex
dynaflex.com
(800) 489-4020

Forestadent USA
forestadentusa.com
(800) 721-4940

Great Lakes Dental Technologies
greatlakesdentaltech.com
(800) 828-7626

Henry Schein Orthodontics
slxclearaligners.com
(800) 837-1552

Orchestrate 3D
orchestrate3d.com
(888) 755-8320

				
EZ-Align	Accusmile	smart moves complete	SLX Clear Aligner System	Orchestrate 3D Clear Aligners
The device is manufactured using a medical grade thermoformed material known as FasTrac.	Forestadent Track Material	Proprietary thermoplastic	SLX Aligners proprietary ClearWear Material provides ultimate clarity and fit to reduce the need for attachments.	PETG Blend
Lab	Both	Lab	Lab	Both
EZ-Align Treatment Review Software (TRS)—proprietary software used with the EZ-Align clear aligner system. Features the latest in 3D viewing technology; presents a variety of tooth movement all the way to the posterior teeth; uses basic geometry for attachments for efficient tooth movement; easy case management; allows doctor to communicate refinements in the treatment plan provided.	Onyx Ceph	Software Interface is a cloud-based platform. Doctors can submit cases, and review, modify, and approve aligner plans for fabrication.	The intuitive Digital Dental Exchange (DDX) SLX Approver Portal streamlines and simplifies the entire case submission and treatment planning process. By providing only relevant patient information to experienced dental-trained technicians, doctors will receive fast and accurate case setups.	Orchestrate 3D Treatment Planning Software (Version 4)
The Systems vary in pricing and each system utilizes the most effective treatment plan with the best value for both the doctor and patient.	Contact company	Tiered pricing starting at \$595	Contact company	Depends on the number of clear aligners
EZ-Align provides a unique pre-aligner expansion option. The EZ-X plus EZ-Align System. The EZ-X is a removable clear lateral development appliance utilizing a nickel titanium, spring loaded memory screw. EZ-X reduces or eliminates the need for enamel reduction and attachments prior to aligner therapy.	Various options available for printing, treatment planning, and data processing to meet the specific needs of your practice.	smart moves aligners are manufactured in, and supported by, U.S.-based technicians. smart moves integrates realistic aligner treatment plans with sophisticated 3D software, and is compatible with any intraoral scanner. Our predictability index provides a baseline for case outcomes so your treatment projections are always realistic.	SLX Clear Aligners are uniquely clear, provide a precision fit, and offer numerous workflow and treatment efficiencies. Henry Schein Orthodontics Attach-Less aligner system incorporates the paradigm-changing SAGITTAL FIRST Philosophy with Carriere MOTION 3D CLEAR Appliances for ultimate aesthetics and efficiency.	Orchestrate 3D is an FDA 510(k) approved orthodontic CAD system that allows the doctor to scan, plan, and 3D print clear aligners and retainers all in-house.


Focus On ALIGNER SYSTEMS

**Rocky Mountain
Orthodontics (RMO)**
orthocaps.com
(800) 525-6375

**SML Space Maintainers
Laboratories**
smglobal.com
(800) 423-3270

**Tortorici Ortho/
Clear Fusion**
clearfusionortho.com
(631) 642-1900

TP Orthodontics Inc
tportho.com
(800) 348-8856



Product name	Orthocaps	Clear Moves	Clear Fusion Hybrid Retainer	Refine Complete Aligner System
What materials is the product made of?	Orthocaps is made of a dual layer polymer with a hard-outer layer and a soft inner layer tray material	Zendura	Clear Thermal reinforced with Orthodontic Acrylic	Thermoplastic polymer
Made in lab or in-house	In-house	Both	Lab	Lab
Software interface	Orthocaps iSetup is a web-based solution that can accommodate both Mac and PCs, and will accept any STL file.	Excocad	N/A	TPRX Portal (http://www.tportho.com/custom-lab-services/)
Orthodontist's cost per case	Varies based on type of case and treatment plan	\$70 per tray; up to 5 trays	\$65 per arch	Treatment times vary; contact company for more info
Special features	Pre-filled attachments, high resolution models which allow for clearer aligner with no striations; ability to choose Rx (MBT, ROTH, etc); treatment evaluations/less refinements = on time finish; ability to choose your arch form; day/night aligners, dual-polymer trays allow for better grip and express tooth movements	Clear Moves are fabricated using accurate computer enhanced tooth movement analysis and special formulated tray materials which provide accurate and predictable tooth movement.	Rigid Clear Retainer	The Refine Complete Aligner System addresses end-of-treatment challenges, such as posterior open bite commonly seen in aligner patients. Treatment options are designed for flexibility, and each case includes your choice of a custom or semi-custom finishing appliance. No advanced training is required—cases can be easily submitted through TPRX.

Information in this guide based on data submitted by product manufacturers. Orthodontic Products strives for accuracy in all data but cannot be held responsible for claims made by manufacturers.

tough on plaque

gentle on orthodontics



only

Platypus[®]

gives you both.

**The ONLY orthodontic
flossing device guaranteed
not to harm braces and
orthodontic hardware!**



Platypus Co.

www.platypusco.com

Request your FREE samples:
samples@platypusco.com
toll free 888.543.6313



In-Ovation® X



In-Ovation® Mini
(Available in Titanium)



In-Ovation® C



In-Ovation® Tubes



In-Ovation® R



No Trace®



In-Ovation® L

In-Ovation runs in our family.

The In-Ovation System is a comprehensive, esthetic and technologically advanced approach to self-ligation. With In-Ovation, you have the ability to choose from a wide range of treatment options to best meet the needs of your practice and your patients, including: reduced chair time and treatment time, fewer and faster adjustments, and the control required for optimal orthodontic results.

Schedule a Free Lunch and Learn!

We'll bring lunch to your office and provide a hands-on demonstration for you and your staff.*
Sign up today! Call 1.800.848.1111 or visit dentsplysirona.com/orthodontics/ortholunch